

User Manual

Mandala Light Alchemy Pulsar




Please read the manual before using the tools.

Version 2.2

Table of Contents

1	General Notes.....	2
1.1	Product Feature and Intend Use	2
1.2	Warning	3
1.3	Warranty	3
1.4	Precautions	4
1.5	Maintenance	4
1.6	Product Includes	4
2	Technical Data	5
2.1	Specifications	5
2.2	Overview and Function.....	6
3	Installation and Operation Instruction	7
3.1	Setup	7
3.2	Switch Status Lights.....	8
3.3	Biofeedback Mode Setup.....	12
3.4	Light Scan Mode Setup	14



1 General Notes

1.1 Product Feature and Intend Use

Thank you for using Light Mandalas Product. Mandala Astral Star have been in good condition before leaving factory.

Mandala Light Alchemy Pulsar is a powerful device which has 16 million solid state light sources programmed to recreate the spectral resonances for libraries of elements in periodic table, trace minerals, amino acids, and outputs the precise colors representing each of them. The new added function – frequency of **Element's Atomic Weight Numbers** allowed the device to have more precise and powerful light colors pulses.

It has useful features such as:

1. Audio mode that can pulse the light according to the **sound frequency** to generate a supportive connection between Light and Sound resonance.
2. Control by Android tablet application with plenty of **Spectroscopic light libraries** and frequency list.
3. Integrated with Mandala Scan to receive biofeedback during meditation.

The frequencies emitted trigger the body to rebalance itself by supplying the nervous, endocrine and immune systems with the information they require. The light frequencies, light nutrients, ex: amino acids, gold, oxygen, etc., assist to bring vitality and well-being back into all the 'bodies' creating a complete balanced being, and clarity toward inner-peace.

Mandala Light Alchemy Pulsar is a complete, all-encompassing 'Pulse, Color, Light and Audio Therapy' tool which makes it excellent for health practitioners, as well as anyone interested in optimal health and emotional wellness. Deep states of meditation can be easily achieved and has been experienced by many users when using **Mandala Light Alchemy Pulsar**.

1.2 Warning

Pay close attention to each safety REMARK and WARNING statement in the manual. This manual contains important operation information.

This product is in perfect condition before delivery, please check all accessory when unpacking the unit.

In order to safely and effectively use this product, please be sure to read this “User Manual” for the proper method usage.

For your eyes safety, please wear light protective goggles and avoid light to expose to your eyes.

1.3 Warranty

Each single product passes through a quality process. Nevertheless, if failures occur, contact the customer service at once. The product is under 12 months (1 year) warranty effectively according to the delivery date stated on the warranty card.

Warranty does not apply to damages, which result from misuse or neglect. The warranty also expires if you dismantle the product or failed to provide the warranty card. The manufacturer is not liable for consequential damage or in case of a non-intended use of the product.

If a failure occurs during the warranty period, the product will be replaced or repaired. The manufacturer reserves the right to exchange components of the product instead of repairing it. If the failure results from misuse or neglect, the user must pay for all cost that may occur. The freight cost and other costs will be depending on which Term & Conditions criteria the product falls into. **Import tax and Customs clearance duty cost is not included in the term of warranty.** In that case, you may ask for an estimated cost beforehand.

1.4 Precautions

1. Please check Mandala Light Alchemy Pulsar for damage during shipping. Do not use the device if any damage is discovered. If the device is damaged, please contact the dealer or Light Mandalas immediately
2. Please Confirm the product you have received are in correct voltage and frequency for your country.
3. Do not disassemble or modify the device in any case.
4. Do not use the device outdoors when rains
5. Places the device in a well-ventilated area.
6. Keep dry and avoid excessive humidity, heat or dust.
7. Do not place the device closer than 2 feet from combustible object.
8. Do not look directly into the tool to avoid harm eyes.
9. Disconnect power when cleaning the tool or not in use.

1.5 Maintenance

1. Clean with soft cloth, humid tissue moistened with water or a water-based glass cleaner.
2. In order to guarantee the fixture works properly, it's essential to keep it clean.

1.6 Product Includes




2 Technical Data

2.1 Specifications

Input Voltage:	AC 110-220 V 50-60 Hz
Power:	100W
Light Source:	12pcs, 10W LEDs (RGB)
LEDs Lifetimes:	100000 hours
Dimmer:	0-100% linear dimmer
communication:	Bluetooth and Audio Bluetooth
Product size:	240 x 230 x 90 mm (W x H x D)
Gross Weight:	1.5 kg

2.2 Overview and Function



3 Installation and Operation Instruction

3.1 Setup

1. Plug in the Power Cord (AC 110-220V, 50/60 Hz only). Make sure the power cord is tightly attached. Then, Turn on the Power Switch.



2. When the power is turned on, the RGB LED light will pulse during the initializing process.



- ***DO NOT look directly to the LED light to avoid harm eyes.***
- ***It's recommended to adjust the LED light facing outward prior to turning ON the device.***

3. After the initializing process, connect the device with Mandala Star App to control the light output.

3.2 Switch Status Lights

There are 3 status lights to select the function for Mandala Light Alchemy Pulsar: Bluetooth mode, Microphone sensitivity mode, Bluetooth audio mode.

1. When the device is turned on, the first light is the red light will always be the Bluetooth light status. This function is needed to connect or pair Bluetooth with Mandala Star Application.



2. Microphone sensitivity mode is the second mode of Light Alchemy, it indicates orange color light. Users can use with Mandala Star Application to insert the healing sound during a healing session. LED light will pulse according to the sound.

Press once to change to the second mode, the status light will be orange color. Then turn the Pulsing Light Brightness Sensitivity Knob to the right until there is almost light coming out of the LED (if the LED light is brighter, turn it slightly back to the left). You can test by clapping your hands, the light will flash in response to clapping. Light Alchemy will flash according to the music being played or follow the sound coming from the speaker.



Push Switch Change mode



Turn the Sensitivity Knob to the right



Clap your hands, Light Alchemy will flash according to your applause.

3. Bluetooth audio mode is the third mode of Light Alchemy, the light Status is Green. LED in this mode, the light can be adjusted and pulse according to audio/music frequency. For Bluetooth Audio, users need to pair Bluetooth from tablet/smartphone setting to connect with another Audio Bluetooth named "XY-ABT".



Push Switch Change mode



LED Status Chang to Green

1. Plug the headset audio cable into the device **Headset Socket**.



2. Use **Mandala Star App** to add your desired music from tablet application and play it. Then users can go back to play the presets from the library. (Please refer to Mandala Astral Star- App User Manual)
3. Adjust the pulsing light brightness with the **Pulsing Light Brightness Sensitivity Knob**.



4. Now you can enjoy Pulsing Color Light therapy as well as Audio/Music therapy at the same time.



3.3 Biofeedback Mode Setup

Biofeedback mode can be used to observe the user's current state during their meditation or spiritual practice. Therefore, users can easily know their current state and allow them to efficiently train their focus.

** Please take note that Biofeedback mode has to be used together with Mandala Dynamic Scan.*

Mandala Dynamic Scan has to be paired with Bluetooth in Mandala Star Application to activate with Light Alchemy.

1. Plug in Ear Clip sensor if user wishes to use this method. Clip it on either one of your earlaps.
2. Use **Mandala Star App** to pair with Light Alchemy and activate the **Biofeedback** pairing from the App setting. (Please refer to Light Alchemy Pulsar- App User Manual)
3. Now the Light Alchemy will project the light according to your current state.

Red: Mental/Emotional level of functioning.



Blue: Heart, mind and soul beginning to connect.



Green: Coherent State. Merge between heart, mind and soul



3.4 Light Scan Mode Setup

In **Light Scan** mode, a suitable light will be projected according to user body need in their current state.

To set up this mode, simply follow the same steps in Biofeedback. Use **Mandala Star App** to pair with Light Alchemy and activate the **Light Scan** mode on Light Alchemy Home Page. (Please refer to Light Alchemy Pulsar- App User Manual)



4. Fault Diagnosis

Failures	Remedies
No power	<ul style="list-style-type: none"> • Please check the power cord if it is tightly attached to the power socket.
No light	<ul style="list-style-type: none"> • Check the tablet Bluetooth setting. Already paired? • Check Mandala Star App setting. Select correct device name?
No pulsing light	<ul style="list-style-type: none"> • Check if the music is playing. • Adjust the Pulsing Light Brightness Sensitivity Knob

Disclaimer

All content found on the Light-mandalas.com Website, Facebook Page, YouTube, including: text, images, video, audio, or other formats are intended for informational, educational and research purposes only and neither Light Mandalas Co. Ltd. nor the publishers of this content take responsibility for possible health consequences of any person or persons reading or following the information in this educational content and is not intended to be a substitute for professional medical advice, diagnosis, treatment or cure. Always seek the advice of your physician or other qualified health professionals with any questions you may have regarding a medical condition.

“No claims other than stress detection and stress reduction may be made”

Light Mandalas products are to be used as a Heart Rate Variability Biofeedback system and Light, Color and Sound therapy, designed for stress detection and stress reduction only. The devices do not evaluate any issue other than stress. Stress can come from many sources; this system uses light and color and sound therapies to reduce stress. Remember, there has never been any reported side effects using light, color and sound therapy anywhere in the world and throughout time.

Notes on the practices of low-level laser therapy (LLLT) and Biofeedback

Low-level laser therapy (LLLT) is considered investigational by the FDA, and no medical claims of cures are permitted in the US; at this time however, low-level lasers are defined within the scope of complimentary medical practices such as physiotherapy, aesthetic, Traditional Chinese Medicine and acupuncture, etc.

The FDA does not regulate Bio resonance tools but they do regulate biofeedback tools and some Bio resonance tools are registered under this category.

*It is a legal requirement at this time that the user does not make any medical claims of diagnosis or cure.

Precaution

Light Mandalas products are working on energetic level. The standard benefit is energy clearance. During the process, your body and mind can be uncomfortable. The symptoms associated with overly rapid clearance of toxins. If those detox symptoms are too strong, you may need to stop using the products. You have to allow the clearance or detox to proceed and you can start to try the products afterwards. However, if the issues still remain the same, it may need to be completely discontinued for your own safety

Shall there be any inquiry, please contact us for further support.

Light Mandalas Co., Ltd. (Thailand Headquarters):

- **Light Mandalas Website:** www.light-mandalas.com
- **Official:** info@light-mandalas.com
- **Instagram:** [@mandalawave9](https://www.instagram.com/mandalawave9)
- **YouTube Channel:** Light Mandalas Official
- **Technical Support Email:** cservice@light-mandalas.com

Local agents contacts:

- **Japan:** japan@light-mandalas.com
- **Taiwan:** lovelight777888@gmail.com
- **Indonesia:** lightseed2020bali@gmail.com
- **West Malaysia:** lightworker@healingtech.org
- **East Malaysia:** malaysia@light-mandalas.com
- **Hungary:** horesnyibela@light-mandalas.com
- **Switzerland, Germany, Austria and the Netherlands:** info@light-mandalas.ch
- **France, Spain:** france-spain@light-mandalas.com
- **USA:** usa@light-mandalas.com
- **New Zealand, Australia:** tachyonhealing.nz@gmail.com